

All leep accessories are designed and built for your particular system AT A SAVINGS OF UP TO 30%



DST Disposable Smoke Tubing

DST (Disposable Smoke Tubes) are designed specifically for LEEP speculum with tubing clips.

#6150K | 50 each | \$50 bag



Stackhouse AirSafe Mini Vac Filter

All-in-one filter. Efficient carbon design removes smoke and odors. Fits Stackhouse, Apen Vac, Valley Vac and Cryomedics systems.

#SFS-300 | 3 each | \$115 box



Patient Tubing Set

Each tube set includes 6' disposable and DST tube. Goes from filter reducer to speculum on Cooper Smoke Evacuation system.

#6084K | 10 each | \$85 box



Stackhouse Smoke Evacuation Tubing

6' smoke evacuation tubing with 24" speculum tube. Fits Stackhouse, Apen Vac, Valley Vac and Cryomedics Airsafe systems.

#SFS-100 | 3 each | \$50 box



Filter Reducer

Goes from the prefilter to patient tubing set on Cooper Smoke Evacuation system.

#6083K | 10 each | \$43 box



Wallach Tube Sets

Replacement set for Wallach Biovac. Includes viral filter, corrugated and speculum tubing.


#9094K | 3 sets | \$85 box



Pre-Filter

Pre-filter for Cooper Smoke Evacuation system.

#6081K | 5 each | \$67.50 box



Patient Return Pads

These single use patient return pads connect from the surgical generator to the patient.

#6050K Single Foil Pad
20 each | \$100 box

#6051K Dual Foil Pad for REM systems
20 each | \$120 box



Uipa Filter

Main filter for the Cooper Smoke Evacuation system for maximum particle capture. Highly efficient carbon design reduces procedural odors. Replace every six months.

#6082K | 1 each | \$110 each



Loop Handpencils

The sterile single patient use handpiece used with any of our disposable electrodes. The wand is operated with a footswitch for ultimate hand precision.

#6041K Finger Controlled
10 each | \$70 box

#6040K Foot Controlled
10 each | \$45 box